

Ordering Information

THD Basic Plus

Code 800075
2 multi-use Doppler probes (10 shots each)
2 reusable needle holders
20 sterile Kits

- Each Kit contains:
- 1 THD device
 - 1 THD Surgy Maxi
 - 1 set of absorbable sutures



THD Basic

Code 800049C
2 multi-use Doppler probes (10 shots each)
2 reusable needle holders
20 sterile Kits

- Each Kit contains:
- 1 THD device
 - 1 THD Surgy Mini



Code 800049D
2 multi-use Doppler probes (10 shots each)
2 reusable needle holders
20 sterile Kits

- Each Kit contains:
- 1 THD device
 - 1 THD Surgy Maxi

Optional



THD Evolution

Code 800000



THD ShiningLight

Code 800051

360118 - 12/2017



THD Basic

The Smart Solution



THD S.p.A.
Via Industria, 1
42015 Correggio (RE) Italy

Customer Service
Ph. +39 0522 634343 - Fax +39 0522 693552
www.thdlab.com



THD Basic

The Smart Solution

DEARTERIALIZATION

Identification and selective ligation of the terminal branches of the superior rectal artery (6 vessels).

Standardised and Reliable

The guide into which the tip of the needle holder is inserted guarantees standardised and repeatable ligation and stitch depth control.



Doppler Guide



Precise ligation



NO excision



THD Basic



THD Basic
CAN BE USED TOGETHER WITH:
THD Surgy Maxi or THD Surgy Mini

Functional Illumination

The design and the optical fibre system provides excellent illumination of the tissue being treated for both THD Basic and THD Surgy.

Optimal Surgical Field

THD Basic and THD Surgy have been designed to offer high functionality and to provide an ideal surgical field and perfect operating access. THD Basic device is optimized for dearterialization, while THD Surgy provides optimal exposure during pexy of the rectal mucosa.

MUCOPEXY

Repositioning of the cushions in their anatomical site by folding the mucohaemorrhoidal prolapse.



THD Surgy Maxi



THD Surgy Mini